Fg-Psd-Pm 05, Version 3, Feb 2022

may not be reused

Ref. PM001-1, PM001-2, PM001-3, PM001-4, PM001-5, PM001-6, PM001-7, PM001-8, PM001-9, PM001-10
Provita Surgical Face Mask Type IIR

Manufacturer/Swedish Office:
 Fantai Group AB
 Örkroken 17
 138 40, Älta, Sweden
 info@provitaswe.se

EU Declaration of Conformity is accessible at:
 www.provitaswe.se

Intended Use

Provita Surgical Face Mask Type IIR is intended to protect against the transmission of infective agents from the users mouth and nose. It is intended to be worn by healthcare personnel during surgical procedures in operating theatres or other medical settings, as well as by patients and members of the general public.

Intended User

The intended users are healthcare personnel, as well as patients and members of the general public.

Indications for use

Provita Surgical Face Mask Type IIR is indicated to be used

- by healthcare personnel when performing, operating or otherwise participating during surgical procedures.
- by patients or persons in the public to reduce the risk of spreading infections, specifically in epidemic or pandemic situations.

Contraindications

Provita Surgical Face Mask Type IIR should not be used by people who are allergic to any material in the device.

Clinical benefits


Provita Surgical Face Mask Type IIR is a source control device that when used can reduce risk of transmission of airborne infections from the user's respiratory tract.

Description

Provita Surgical Face Mask Type IIR is made of non-woven polypropylene and material in 3-layers with elastic loop for ears. The

face mask is 100% latex free. >98% bacterial filtration efficiency. Dispose according to hospital policy and local regulation.

The face masks are made of following material:

1. Nose clip material: Polypropylene and steel.
2. Outer layer: Polypropylene.
3. Second layer: Meltblown fabric.
4. Inner material: Polypropylene.
5. Ear loops: Polyurethane.

Performance

This product meets the requirements of the European Standard EN 14683:2019+AC:2019: Surgical face masks, type IIR.

Test

| | |
|--|--------|
| Bacterial filtration efficiency (BFE), (%) | ≥ 98 |
| Differential pressure (Pa/cm) | < 60 |
| Splash resistance pressure (kPa) | ≥ 16.0 |
| Microbial cleanliness (cfu/g) | ≤ 30 |

Material Performance

The performance of the face masks quoted above is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, mechanical wear etc.

Storage

The face masks must be stored dry and under clean conditions where the temperature is between 10°C–35°C with humidity between 45%–65%. Preferably stored in original packaging under dry conditions.

When stored as recommended, the chemical properties on the non-woven material and the ear loops will remain unchanged for up to 2 years from the date of manufacture.


Warnings and Precautions

- Components used in manufacturing may cause allergic reactions in some users.
- If any allergic reaction occurs, seek medical advice or visit a doctor immediately.
- For single use only. Do not re-use.
- Before usage, inspect face mask for any defect or imperfections.
- Do not use if packaging or face mask is damaged.
- Do not touch the inside of the mask.
- Do not use after expiry date.
- Children or other people in need of assistance must be assisted when applying the face mask.
- Replace the face mask immediately if it is difficult to breathe, it is damaged or fails to maintain a tight fit.

User guide

1 Always remember hand hygiene! Clean your hands before removing the face mask.


2 Make sure the face mask is intact.

3 Take the elastic band and place the face mask over the nose and mouth with text on the outside.

4 Attach the elastic band behind the ears.

5 Pull the upper edge over the back of the nose, and the lower edge under the jaw.


6 Press the bridge of the nose so that the face mask is close to the skin.


7 Make sure the edges of the face mask are close to the skin throughout the circumference.

8 The face mask can be used up to 4 hours.

9 Remove the mask by holding the ear loops.

10 Throw it in a dustbin with a lid, wash your hands.

In case of incident

Any serious incident that occurs in relation to the face mask should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. The manufacturer cannot be held responsible for improper use of the face mask.

 MD stands for Medical Device.

 CE Mark conforms to essential requirements of the MD Regulation.

 Batch code indicates the manufacturer's batch code so that the batch or lot can be identified.

 Keep dry, protect from moisture.

 Temperature range.

 Humidity range.


 MD intended for one use on a single patient during a single procedure.

 Keep away from sunlight.

 Use-by date, the date after which the medical device is not to be used.

 Date of manufacture, indicates the date when the MD was manufactured.

 Manufacturer, indicates the MD manufacturer, as defined in MD Directive.

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