



Respirator BM-027, BM-027V. Protection-class FFP2 NR D is made for environments in which harmful particles can be found in the air we breathe. Respirator masks of this class provide a barrier to at least 94% of particles measuring larger than 0.3 micron and may be used in environments where hazardous substances exceed the OEL by up to a maximum of tenfold concentration. Our filter system, with layer of meltblown cloth, protects wearers from unpleasant odours and particles on top of the required breathing protection.



***Extra optional EAR hook stripe system**



Standard: EN149:2001+A1:2009, FFP2

Options: **BM-027** Strap behind ears. Behind the wearer's head* (*Includes Head Strip Holder)

BM-027V With valve. Strap behind ears. Behind the wearer's head* (*Includes Head Strip Holder)

Filtration efficiency: ≥ 94%

Designation: NR - Non reusable (single shift use only)

Marking: PPE marking, FFP2 NR D, CE1463, EN 149:2001+A1:2009

Materials:

50gsm Spunbond Nonwoven

40gsm Meltblown Nonwoven

40gsm Hot air Nonwoven


25gsm Meltblown Nonwoven

Mask weight: **BM-027** - 5 grams

BM-027V - 7.5 grams


Shelf life: Two years

Packaging BM-027




2 pcs in a polybag, 10 pcs in a inner box, 48 inner boxes in a carton box

- 1 **Size:** 150mm (y) x 64mm (x) x 120mm (z)
Quantity: 1box (10pcs.)
- 2 **Size:** 590mm (y) x 380mm (x) x 280mm (z)
Quantity: 1/48box (480pcs.)
- 3 **Size:** 1180mm (y) x 760mm (x) x 1680mm (z)
Quantity: 24/48box(11.520pcs.)



2 pcs in a polybag, 270 polybags in a carton box.

- 1 **Size:** 220mm (y) x 110mm (x) x 7mm (z)
Quantity: 1pack (2pcs.)
- 2 **Size:** 590mm (y) x 380mm (x) x 280mm (z)
Quantity: 1/270pack (540pcs.)
- 3 **Size:** 1180mm (y) x 760mm (x) x 1680mm (z)
Quantity: 24/270pack (12.960pcs.)




10 pcs in a polybag, 80 polybags in a carton box.

- 1 **Size:** 185mm (y) x 240mm (x) x 25mm (z)
Quantity: 1pack (10pcs.)
- 2 **Size:** 590mm (y) x 380mm (x) x 280mm (z)
Quantity: 1/80pack (800pcs.)
- 3 **Size:** 1180mm (y) x 760mm (x) x 1680mm (z)
Quantity: 24/80pack (19.200pcs.)

Storage Conditions and Shelf Life The filter masks have to be stored in their packets to retain their properties. Shelf life of unopened product is 2 years from date of manufacture when stored within temperature range of -10°C to +40°C and at less than 85% relative humidity.


Packaging BM-027V

1 pcs in a polybag, 10 pcs in a inner box, 36 inner boxes in a carton box




- 1 **Size:** 145mm (y) x 180mm (x) X 126mm (z)
Quantity: 1box (10pcs.)
- 2 **Size:** 590mm (y) x 380mm (x) x 280mm(z)
Quantity: 1/36box (360pcs.)
- 3 **Size:** 1180mm (y) x 760mm (x) x 1680mm(z)
Quantity: 24/360pack (8.640pcs.)

2 pcs in a polybag, 270 polybags in a carton box.



- 1 **Size:** 220mm (y) x 110mm (x) x 7mm (z)
Quantity: 1pack (2pcs.)
- 2 **Size:** 590mm (y) x 380mm (x) x 280mm (z)
Quantity: 1/270pack (540pcs.)
- 3 **Size:** 1180mm (y) x 760mm (x) x 1680mm (z)
Quantity: 24/270pack (12.960pcs.)

Storage Conditions and Shelf Life The filter masks have to be stored in their packets to retain their properties. Shelf life of unopened product is 2 years from date of manufacture when stored within temperature range of -10°C to +40°C and at less than 85% relative humidity.

ENG**BM-027 BM-027V****FFP2 NR D** **Usage instruction**

 Do not wear with beards or other facial hair that may inhibit contact between the face and edge of the respirator


Donning instructions to be followed each time you use the mask. Before use, check for visible damage, damaged or dirty (on breathing side) filtering half masks cannot be used.

1. Open the mask with both hands to the two sides.
2. Straighten the nose bridge stripe.
3. Fit the mask on your face covering the mouth and nose. The straps should not be between the mask and your face.
4. Pull the straps above your ears. Adjust the mask until it fits comfortably.
5. **Adjust the nose bridge stripe, until it completely fits to your nose bridge and you don't feel any air leakage.**
6. **Cover the mask with both hands until the mask fits tightly around the face.**

Attention!!! The fit is critical to the performance of the product!!!

 **Precaution**

1. This product is not suitable for protection against harmful and oily aerosols.
2. The product does not generate oxygen and is not suitable for use in anoxic environments, underwater operations, during evacuation or firefighting environments.
3. It is not recommended to use this protective mask if you do not feel a good adhesion between the face and the edge of the mask during use.
4. If found during use that the mask is damaged or the breathing resistance is too large, replace with a new one.
5. This product is for one-time use, masks should not be reused. Recommended use time is 4-8 hours.
6. It is unlikely that the requirements for leakage will be achieved if facial hair passes under the face seal.

 **Warning - this product is not flame resistant and must not be used in areas with open flames**

Materials used: 4 layer -50gsm Spunbond Nonwoven, 40gsm Meltblown Nonwoven, 40gsm Hot air Nonwoven, 25gsm Meltblown Nonwoven

Storage:    Temperature range of storage conditions Humidity range of storage conditions Protect from direct contact with water	Types of packaging: 1 / 2 / 5 / 10 pcs.
---	---

 BALTIC MASKS, UAB Savanoriu ave. 178F, Vilnius, Lithuania

CE notified body: Polski Rejestr Statków S.A(NB ID: 1463)

These instructions for use have been prepared in accordance with the requirements set out in the Regulation (EU) 2016/425 of the European Parliament and of the Council. [Link to the declaration of conformity and CE certificate: www.balticmasks.com](#)