



Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**



Product description

Order number	62.548.004
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap assembled, with print, label/print: blue/white, with graduations, free from DNA, DNase/RNase, pyrogens/endotoxins, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Type of cap	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Size

Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	free from DNA, DNase/RNase, pyrogens/endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of inner packaging	bag
Piece(s) / inner box	25
Piece(s) / case	300
Piece(s) / pallet	9000
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0128 kg
Weight of case	4.61 kg
EAN of inner box	4038917265201
EAN case	4038917265157

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

