



Product description

Order number	72.380
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, HD-PE, white, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Type of cap	screw cap
Internally or externally threaded cap	Internal thread
Design	with silicon O-ring
Volume of work	1,6 ml
Tube	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	1.6 ml
Max. volume	2 ml
Diameter	12.65 mm
Length of product	47.7 mm
Length excluding cap	37.8 mm

Material & colours

Cap material	High Density Polyethylene (HD-PE)
Colour of cap	white
O-ring material	silicone
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0019 kg
Weight of case	4.93 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
CryoPure tubes, 2 ml, QuickSeal screw cap, white

Page 3

EAN of inner box 4038917081405
EAN case 4038917080606

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT