




**Product Specification**  
**Faeces tube, screw cap, (LxØ): 101 x 16.5 mm,**  
**transparent, sterile**

**Page 1**



**Product description**

Order number	80.623
Product description	Faeces tube, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, sterile

**Product characteristics**

Type of cap	screw cap
-------------	-----------

**Size**

Diameter	16.5 mm
Length of product	102 mm
Length excluding cap	101 mm

**Material & colours**

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Faeces tube, screw cap, (LxØ): 101 x 16.5 mm,**  
**transparent, sterile**

**Page 2**

**Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0068 kg
Weight of case	3.89 kg
EAN of inner box	4038917082983
EAN case	4038917001960

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**