



**Product Specification**  
**Faeces tube, screw cap, (LxØ): 76 x 20 mm,**  
**transparent**

**Page 1**



### Product description

Order number	80.734
Product description	Faeces tube, transparent, screw cap, cap: brown, (LxØ) without cap: 76 x 20 mm

### Product characteristics

Type of cap	screw cap
-------------	-----------

### Size

Diameter	20 mm
Length of product	77 mm
Length excluding cap	76 mm

### Material & colours

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Faeces tube, screw cap, (LxØ): 76 x 20 mm,**  
**transparent**

**Page 2**

**Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	21000
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.007 kg
Weight of case	4.07 kg
EAN of inner box	4038917083812
EAN case	4038917002707

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**