



Product description

Order number	86.1562.010
Product description	Inoculation loop, volume: 10 µl, material: PS, blue, sterile, 10 piece(s)/blister

Size

Width of product	6 mm
Length of product	200 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	blue

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1000
Type of inner packaging	blister
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	132000
Depth of box	216 mm
Width of box	161 mm
Height of box	101 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.0012 kg
Weight of case	1.79 kg
EAN of inner box	4038917130875
EAN case	4038917023337

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

