



Product description

Order number 86.1569.005
Product description Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 4 piece(s)/blister

Product characteristics

Example of use for spreading a bacterial suspension

Size

Width of product 40 mm
Length of product 185 mm

Material & colours

Product material Polystyrene (PS)
Colour of product blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	blister
Piece(s) / inner box	4
Piece(s) / case	500
Piece(s) / pallet	33000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.002 kg
Weight of case	1.7 kg
EAN of inner box	4038917130998
EAN case	4038917023566

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

