



Product description

Order number	20.1331.100
Product description	Microvette® APT 250 K2E, capillary blood collection, preparation: K2 EDTA, K2E, nominal volume: 250 µl, Ø: 13 mm, with plastic label, cap violet, cap, with pierceable membrane, round base, inner receptacle: cylindrical, incl. quick-release cap, non-sterile, 50 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	K2 EDTA
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	transparent/violet
Cap	membrane cap
Type of cap	cap
Application	capillary blood collection
Base shape	round base
Volume of work	250 µl - 250 µl
Shape of inner tube	cylindrical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	250 µl
Diameter	13 mm
Length excluding cap	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Polyethylene (PE)
Colour of cap	violet

Purity & certification

Satisfies the requirement	P 650, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	32000
Depth of case	378 mm
Width of case	268 mm
Height of case	221 mm
Case volume	0.0224 cbm
Weight of product	0.0075 kg
Weight of case	4.01 kg
EAN of inner box	4038917543026
EAN case	4038917543019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

