



**Product Specification**  
**Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP**

**Page 1**



### Product description

Order number	75.9922.744
Product description	Multi-purpose container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

### Product characteristics

Graduation	yes
Type of cap	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	70 ml

### Size

Sample volume	70 ml
Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP**

**Page 2**

Length of product 55 mm  
Length excluding cap 55 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Cap material High Density Polyethylene (HD-PE)  
Colour of cap yellow

**Purity & certification**

Product category In-vitro diagnostic, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of inner packaging bag  
Piece(s) / inner box 250  
Piece(s) / case 500  
Piece(s) / pallet 8000  
Case volume 0.0889 cbm  
Weight of product 0.0123 kg  
Weight of case 7 kg  
EAN of inner box 4038917152754  
EAN case 4038917027793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**