



Product description

Order number	72.690.001
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	39 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Polypropylene (PP)
Colour of cap	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.001 kg
Weight of case	5.56 kg
EAN of inner box	4038917161657
EAN case	4038917554138

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

