



Product description

Order number	60.597.001
Product description	Screw cap tube, working volume: 120 ml, (LxØ): 114 x 44 mm, material: PP, flat base, transparent, PP, cap: natural, cap assembled, sterile, 25 piece(s)/bag

Product characteristics

Graduation	no
Type of cap	screw cap
Base shape	flat base
Volume of work	120 ml

Size

Diameter	44 mm
Length of product	114 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Polypropylene (PP)
Colour of cap	natural

Purity & certification

Satisfies the requirement	ADR
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	100
Type of inner packaging	bag
Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	3600
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0222 kg
Weight of case	2.76 kg
EAN of inner box	4038917114172
EAN case	4038917014212

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

